

QUICK REFERENCE INSTRUCTIONS

Step 1



Put a small amount of substance into a container. Add a small amount of water (a teaspoon) and mix.

The amount should cover the president's hair on a penny.



Step 2




Dip the end of the strip for 15 seconds.

Lay the strip flat.

Step 3

Read results at 5 minutes.



Interpretation of results

NEGATIVE	POSITIVE	INVALID
(No Benzodiazepines Detected)	(Benzodiazepines Detected)	(Perform Test Again)

WARNING: The results of the SAFe-B test are not a guarantee that the tested substance can be safely consumed. Even if the results of the SAFe-B test are negative for benzodiazepines, **anyone consuming the tested substance could overdose, sustain severe injury, or die.** Failure to follow the above instructions could result in false positive or false negative test results. The above instructions are only for the SAFe-B test, and you should not use them for any tests other than the SAFe-B test. **Please see the reverse side for other warnings and information regarding the SAFe-B test.**

Please see other side for Quick Reference Instructions

Warnings

- Read, understand, and follow the instructions on the reverse side, your failure to do so could result in false positive or false negative test results, severe injury or death.
- Even if the results of the SAFe-B test are negative for benzodiazepines, anyone consuming the tested substance could overdose, sustain severe injury, or die.
- Do not consume illicit drugs, especially if you suspect that they contain, or the SAFe-B test result is positive for, benzodiazepines.
- The SAFe-B test does not detect the presence of substances other than benzodiazepines, and it does not detect all forms of benzodiazepines.
- The SAFe-B test does not detect the percentage of benzodiazepines present in the tested substance.
- Exposure to even small amounts of benzodiazepines, or handling of benzodiazepines, can cause severe injury or death.
- When testing a substance, use protective equipment, such as gloves and eye protection.
- Do not use the SAFe-B test after expiration or if the packaging is damaged or shows other signs of tampering.
- Dispose of the SAFe-B test in accordance with applicable law.
- Anyone using the SAFe-B test accepts all responsibility for any injury, death, criminal penalty, or civil liability that occurs from consumption of a substance, whether an illicit drug or not, and regardless of whether the substance was tested, or was not tested, for benzodiazepines.
- Call 911 or seek immediate medical attention if you think you or someone else has overdosed or is at risk of overdose. Common signs of an overdose are listed in the "How do I recognize an overdose?" section below.
- If you or someone you know requires counseling or treatment for drug use, call the SAMSHA National Helpline at 1(800)-622-HELP(4357). This helpline is confidential and available 24/7, 365 days a year.

How do I recognize an overdose?

Signs of an overdose include:

- Small, constricted "pinpoint pupils"
- Falling asleep or loss of consciousness. Slow, shallow breathing
- Choking or gurgling sounds
- Limp body
- Pale, blue, or cold skin
- Because this list is not exhaustive, call 911 or seek immediate medication attention if you think you or someone else has overdosed or is at risk of overdose.