

# The Difficult Airway

## Strategy and Tactics

Don Byars, MD, FACEP  
Assistant Professor of Emergency Medicine  
Eastern Virginia Medical School

Operational Medical Director  
City of Portsmouth Department of Fire and EMS  
City of Portsmouth Police Department TRU

# Objectives

- Identify the Enemy
- Outline Our Strategy
- Define Our Tactical Options
- Review and Demonstrate New Tactics
  - KING Airway
  - I-LMA
  - Video Laryngoscopes

# The Difficult Airway

- Rule #1
- ALL airways must be treated as difficult airways!




# A Few Definitions

- *Difficult Airway*
- *Crash Airway*
- *Failed Airway*

# Identify the Enemy (Airway)


- Various Assessment Schemes
  - 3-3-2
  - Mallampati Score
  - Facial Structure
  - Tracheotomy Scars
  - Facial Hair
  - Etc, etc, etc



# Rapid Fire Assessment

- 10 Second Difficult Airway Assessment
  - “LEMON” Approach
    - L-ook for obvious signs of potential problems
    - E-valueate via the 3-3-2 rule
    - M-allampati
    - O-bstruction
    - N-eck Mobility

# 3-3-2 Rule




# Mallampati




Class I: soft palate, uvula, fauces, pillars visible  
No difficulty



Class II: soft palate, uvula, fauces visible  
No difficulty



Class III: soft palate, base of uvula visible  
Moderate difficulty



Class IV: hard palate only visible  
Severe difficulty

# Outline Our Strategy

- Strategy:
- An Overall plan of action designed to accomplish a specific goal



# The Oft Forgotten...

## BVM

- More lives are saved from solid bag valve mask technique than from any so-called “advanced airway device” ever made
- All airway classes should begin and end with BVM technique


# Approach to the Difficult Airway

- Traditional “Double-Setup”
  1. Laryngoscope and tube
  2. Cricothyrotomy Tray (kit)



# Approach to the Difficult Airway

- “Deca-Setup”
  1. Traditional Laryngoscope
  2. GEB
  3. I-LMA
  4. Glidescope
  5. Cricothyrotomy Kit
  6. KING
  7. Retro-grade Kit
  8. Fiber Optic Scope
  9. Trans-Tracheal Jet Ventilator
  10. Tracheostomy Tray




# “Awake” Technique?

- Purpose: to determine whether direct laryngoscopy is **possible!**
- What does awake really mean?
- The patient maintains spontaneous respirations
- Sedation + Topical Anesthesia

# “Awake” Technique

- This is the classic case for an “awake” nasal fiberoptic intubation
- KEY point is to NOT paralyze in the field!



# Grade I View

- Grade I View
- Whole Larynx Visible



# Grade 2 View

- Grade 2 View
- Only Part of the Glottis Visible



# Grade 3 View

- Grade 3 View
- Epiglottis Only



# Grade 4 View

- Grade 4 View
- NO laryngeal structures




# Awake Technique

- Grade 1 view = RSI ok
- Grade 2 view = RSI ok
- Grade 3 view = NO RSI
- Grade 4 view = “Impossible” Intubation

# Definition

- Failed Airway
- 3 Failed Attempts by an experienced operator using optimal patient positioning and best possible technique
- O2 Sat < 90% with 2 man BVM + adjunct
- Grade IV view on Laryngoscopy

- Failed Airway Algorithm



# Define Our Tactical Options

- Tactics:
- Near term actions taken to solve a specific problem or achieve a certain aim




# Difficult Airway Options

- GEB
- KING
- I-LMA
- Glidescope
- Cricothyrotomy Kit
- ~~● Fiberoptic Scope~~

# Gum Elastic Bougie


- Tracheal Tube Introducer
- 60 cm long
- 15 french diameter
- J-tip curve at tip



# Set Up

- Set up for traditional intubation as usual
  - Same equipment
  - Same procedure
  - Same positioning
  - Same drugs
- Prepare GEB
  - J-tip
  - Gentle C curve






# Attempt Intubation

- Consider GEB IF:
  - Grade II or III view seen (Arytenoids or epiglottis only)
- Obtain best possible view with laryngoscope



# Insert GEB

- If cords **NOT** visible pass GEB so that J-tip points anterior when passing the epiglottis



# Is It In?

- Characteristic “clicks” should be felt when passing the GEB over tracheal rings
- If unsure, advance the GEB to see if resistance is met at the carina
- **IF CLICKS ARE NOT FELT AND NO RESISTANCE IS MET WITH ADVANCEMENT OF THE GEB => ESOPHAGEAL PLACEMENT**

# Confirmation

- Once the tube is passed through the cords, the intubator will maintain control of the ETT
- The assistant will remove the GEB
- Endotracheal tube placement will be confirmed using traditional methods

# Questions on the GEB?



# KING LTD Airway

- THE prehospital rescue airway device
- BLS skill level for insertion technique
- 10 seconds to oxygenation



# KING LTD Airway

- Two Balloons
- One in the distal esophagus
- Other in the posterior pharynx
- Ventilation outlets



# **KING LT Insertion**

# KING LTD Airway

- Chin Lift
- Tube held “cocked” at 45 degree angle
- CAREFUL with the balloons going past the teeth!!!



# KING LTD Airway

- Rotate the tube back to neutral as the tube tip passes the tongue



# KING LTD Airway

- Advance the tube **WITHOUT** force



# KING LTD Airway


- Until the tube phalange is resting on the teeth or gums
- Then inflate cuff with 60 cc of air
  - or to 60 cm H<sub>2</sub>O



# Questions on the KING LTD?



# The LMA Fastrach



# The LMA Fastrach

- Facilitates oxygenation and ventilation with or without intubation
- LMA portion alone is NOT a definitive airway
- Accommodates up to an 8.0 ET tube
- Inserted blindly

# Tracheal Intubation Using The Fastrach

© 2000/2001

---


---

# Insertion Technique


- Visually inspect to ensure device is not defective/damaged
- Deflate cuff fully – will not hold tight with leaky cuff
- Place patient's head in neutral position
- Lubricate POSTERIOR tip

# Insertion technique

- Place tip against hard palate with metal tube against chin




- Using a circular motion, glide the mask along the hard palate and into position



# Insertion Technique

- Inflate cuff using 30 cc syringe to approx 50% of maximum
  - 15 to 20 cc of air
- Use fine adjustments of mask position to find “best fit” with minimal air leak



# Insertion Technique

- Use end-tidal CO<sub>2</sub> detector to confirm placement
- Oxygenate with bag-valve mask 15 L/  
Min oxygen
- Intubate




**FIGURE 8**


*Ventilate the patient prior to intubation.*

# Intubating through the I-LMA

- Inspect the ET tube
- Lubricate ET tube
- Use tube to lubricate LMA



# Intubating through the I-LMA



# Intubating through the LMA

- Lift straight up on handle to ensure proper LMA seal
  - The “Chandy” maneuver
- Gently pass tube to minimum depth of 15 cm



# Intubating through the LMA

- Inflate ETT cuff using 10 cc syringe
- Confirm placement with end-tidal CO<sub>2</sub> detector, auscultation for breath sounds, watch O<sub>2</sub> sats





# Questions on the I-LMA?




# Cricothyrotomy

- The ONLY option for a true FAILED airway
- Can't Intubate + Cant Ventilate = Knife
- Multiple Options to accomplish the goal

# Cricothyrotomy



# Cricothyrotomy



**The Bougie-Aided Cricothyrotomy**  
with Darren Braude, MD, EMT-P

# Questions on the GEB Cric?



# Glidescope

- Video Laryngoscopes ARE THE FUTURE!!
- Provide Grade I view 99% of the time
  - In trained hands



# Glidescope Ranger

- 7" by 6"
- Weighs 1.5 lbs
- Can survive extreme temperatures



# Glidescope

- Field tested in Iraq
- Used primarily at Combat Support Hospitals



# Glidescope



434

445

474

"Rescue"

**KING**

**VS**

**Video Laryngoscope**

THE END

