
Clinic Operations Scheduled Duties Mobile Units

Daily

Turn on all systems and check
Bag instrument set ups
Clean, lubricate and sterilize handpieces
Clean, disinfect, sterilize instrumentation
Flush unit hoses
Rinse suction system with Purevac cleaner or equivalent product
Fill unit water reservoir with fresh water and add ICX or equivalent

Weekly

Restock all supplies
Change ultrasonic solution
Clean water line hoses with Sterilex Ultra or equivalent per manufacturer's instructions
Spore test autoclave and document temperature and pressure

Monthly

Complete supply orders
Clean autoclave as needed
Clean or replace unit filters

Annually

Test delivery unit water quality and document

10/02/2018