

Laboratory Testing for Pertussis

Test	Sensitivity	Specificity	Comments
DFA	30-71%	66-100%	Should be used only as a screening test in conjunction with culture and PCR
Culture	11-73%	100%	Organism is extremely labile; requires prompt transport and plating; sensitivity best when specimen collected during catarrhal stage of the disease
PCR	61-99%	88-100%	No single technique is universally accepted or validated; no FDA cleared test available; culture should also be performed whenever PCR testing is performed
Serological Tests			Difficult to interpret, especially in immunized people; no FDA cleared test available; lack of standardized interpretive criteria
IgM			Generally not helpful
IgG	63%*	99%*	Specificity depends on the antigen used (PT is specific); more sensitive than IgA; paired sera can also be tested
IgA			Specificity depends on the antigen used; rarely detected in uninfected vaccinees; young infants may not be able to make IgA

References:

- Centers for Disease Control and Prevention. 2000. Guidelines for the Control of Pertussis Outbreaks. CDC Atlanta, GA.
- Chan EL, et al. 2002. The Use of TaqMan PCR Assay for Detection of *Bordetella pertussis* Infection from Clinical Specimens. Arch Pathol Lab Med. 126: 173-176
- Loeffelholz MJ, Thompson CJ, Lon, KS, Gilchrist MJR. 1999. Comparison of PCR, Culture, and Direct Fluorescent-Antibody Testing for Detection of *Bordetella pertussis*. J Clin Micro 37: 2872-2876
- McGowan KL. 2002. Diagnostic Tests for Pertussis: Culture vs. DFA vs. PCR. Clin Micro Newsletter 24: 143-150
- Muller F-MC, Hoppe JE, vonKonig C-HW. 1997. Minireview: Laboratory Diagnosis of Pertussis: State of the Art in 1997. J Clin Micro 35: 2435-2443

*Massachusetts State Laboratory Institute's test, the only serological test recognized by CDC as confirmatory and only for patients in Massachusetts